

STAY SAFE with our PROTECTION PACK

Face Shield

For Professional Protection

FEATURES & BENEFITS



Protective from

- Premium foam for forehead provides extra comfort
- · Wrap-round design for full protection
- High temperature and shock resistance
- Excellent anti-fog performance on both side
- Easy wearing
- Comfort and soft forehead protection
- Tight and stable holding

MEASUREMENT & MATERIAL

- 210(h) x 160(w)
- 100% Transparent PVC
- No burr & sharp edges





FDA IR Thermometer

Designed to scan and monitor body temperature









SPECIFICATIONS

- Effective Distance: 5-8cm
- Measurement Method: Non-contact
- Memory: Storage of last 10 readings
- Display: 3 Colors Backlight LCD (red, green, orange)

tGuard

Fever Screening System

tGuard adopts an infrared radiation measurement technique to measure body temperature by non-contact means. It provides audio alarms when a person with elevated temperature passes by.

KEY SPECIFICATION



Measurement Distance 0.15 - 2m



High Accuracy ±0.5°C ±0.3°C (with blackbody)



Wafer Level Packaging



Temperature Measurement Range



Recording and Searching



Auto Audio Alarm





TG-2200



YEAR

COMPONENTS











TG-2300

Tablet Power Line

Mini-USB to USB-A

USB-A to USB-C

Infrared Camera



Extended Adapting Piece



Metallic Support for Tablet



V

NOVATEC SANITARYWARE LANKA (PVT) LTDNo. 9, Chilaw Road, Negombo,
Sri Lanka



+94 31 4936843/31 4936847



sales.lk@sanieuro.com sales.novateclanka@gmail.com



www.sanieuro.com/lk