

STAY SAFE with our PROTECTION PACK

Face Shield

For Professional Protection

FEATURES & BENEFITS



Protective from

- Premium foam for forehead provides extra comfort
- · Wrap-round design for full protection
- High temperature and shock resistance
- Excellent anti-fog performance on both side
- Easy wearing
- Comfort and soft forehead protection
- Tight and stable holding

MEASUREMENT & MATERIAL

- 210(h) x 160(w)
- 100% Transparent PVC
- No burr & sharp edges





FDA IR Thermometer

Designed to scan and monitor body temperature









SPECIFICATIONS

- Effective Distance: 5-8cm
- Measurement Method: Non-contact
- Memory: Storage of last 10 readings
- Display: 3 Colors Backlight LCD (red, green, orange)

tGuard

Fever Screening System

tGuard adopts an infrared radiation measurement technique to measure body temperature by non-contact means. It provides audio alarms when a person with elevated temperature passes by.

KEY SPECIFICATION



Measurement Distance 0.15 - 2m



High Accuracy ±0.5°C ±0.3°C (with blackbody)



Wafer Level **Packaging**



Temperature Measurement Range



Recording and Searching



Auto Audio





TG-2200





TG-2300

tGuard Lite

COMPONENTS











Mini-USB to USB-A

USB-A to USB-C

Metallic Support

for Tablet

Infrared Camera



Tablet Power Line

Extended Adapting Piece

SANIEURO MYANMAR CO., LTD.

No. 12, Room 7, Parami Road, Chaw Dwin Gone, Ma Yan Gone Township, 11061 Yangon, Myanmar







sales.mm@sanieuro.com



www.sanieuro.com/mm